
Isolation of drug target protein (3) Histone Deacetylase Inhibitor【Vorinostat】

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

Fig.1 Immobilization of Vorinostat amino derivative on NHS beads

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Materials and method

Materials

1. Vorinostat amino derivative
2. NHS beads
3. HeLa cell extracts (cytosolic fraction) – 3mg/ml
4. Binding & Washing Buffer (20mM HEPES-NaOH(pH7.9), 100mM KCl, 1mM MgCl₂, 0.2mM CaCl₂, 0.2mM EDTA, 10%(v/v) glycerol, 0.1% NP-40, 1mM DTT, 0.2mM PMSF)
5. Elution Buffer (0.0625M Tris-HCl (pH6.8), 0.005% BPB, 2% SDS, 10% glycerol, 5% 2-mercapto ethanol)
6. Vorinostat
7. Rabbit polyclonal anti-HDAC1 (from Santa Cruz Biotechnology)
8. Rabbit polyclonal anti-HDAC2 (from Santa Cruz Biotechnology)
9. Rabbit polyclonal anti-HDAC3 (from Santa Cruz Biotechnology)
10. Rabbit polyclonal anti-HDAC8 (from Santa Cruz Biotechnology)
11. Anti-Rabbit IgG, HRP-Linked Whole Ab Sheep (from GE Healthcare)
12. Starting Block™ (TBS) Blocking Buffer (from Thermo scientific)
13. TBS-T buffer (137mM NaCl, 2.68mM KCl, 25mM Tris-HCl, pH7.4, 0.1w/v% Tween-20)

Methods 1 (Immobilization)

1. Apply
 - 0.5mg NHS beads and 0.3mM vorinostat amino derivative in 100ul DMF
 - 0.5mg NHS beads and NO vorinostat amino derivative in 100ul DMF (Negative Control)
2. Reaction
 - 1) Immobilization
 - 70min at r.t.
 - 2) Masking
 - Inactivate unreacted NHS according to recommended conditions.

Methods 2 (Affinity purification)

1. Wash
 - Wash beads with washing buffer 3 times at 4°C (on ice).
2. Add sample solution
 - Add 1000ul HeLa cell extracts to beads.
3. Reaction
 - Incubate for 240min at 4°C.
4. Wash
 - Remove sample solution.
 - Wash beads with washing buffer 3 times at 4°C (on ice).
5. Elution
 - Vorinostat elution : -
 - Add 40ul elution buffer and resuspend beads.
 - Boil for 5min at 98°C and remove beads.
 - Vorinostat elution : +
 - Add 30ul of 25mM vorinostat in washing buffer and resuspend beads.
 - Place the tube on ice for an hour to allow bound proteins to elute by tapping at intervals.
 - After spin down, separate magnetically.
 - Transfer the supernatant to a new clean tube.
 - Add 10ul of 4 × elution buffer to the tube.
 - Boil for 5min at 98°C.
6. Analyze the samples by SDS-PAGE and silver staining

FG beads information

Product name	NHS beads
Product number	TAS8848N1141
Storage temperature	-20°C
Storage buffer	Isopropyl alcohol
Size of beads	190nm ± 20nm
Functional groups	200 - 300nmol/mg

Methods 3 (Western blotting)

1. Perform SDS-PAGE and transfer the protein to a PVDF membrane.
2. Block the membrane with the blocking buffer for 15min at room temperature.
3. Dilute the anti-HDAC antibody with the blocking buffer to 1/200.
4. Incubate for 60min at room temperature.
5. Dilute the secondary antibody with the TBS-T buffer to 1/1000.
6. Incubate for 60min at room temperature.
7. Wash the membrane with TBS-T buffer three times.
8. Detect with a chemiluminescence substrate.

Methods 4 (MS)

1. Excising the bands from the gel.
2. Wash the band pieces, then dry the band pieces.
3. Add 500ul/band of 10mM DTT/25mM NH₄HCO₃.
4. Incubate for 60min at 56°C.
5. Add 500ul/band of 55mM iodeacetamide/25mM NH₄HCO₃.
6. Incubate for 60min at room temperature.
7. Wash the band pieces with 200ul of 25mM NH₄HCO₃ for 10min.
8. Wash the band pieces with 500ul of 25mM CH₃CN/NH₄HCO₃ (50:50 v/v) for 10min.
9. Dry the band pieces.
10. Add 10ul/band of 10ug/ml trypsin in 25mM NH₄HCO₃.
11. Incubate overnight at 37°C.
12. Add 50ul/band of 25mM CH₃CN/trifluoroacetic acid/ultrapure water (50:5:45 v/v).
13. Incubate for 30min at room temperature.
14. Transfer supernatant to a new clean tube.
15. Dry the peptides in the tube.
16. Injection on the MS.