

Experiment Protocol 401

Preparation of cell extract (Large scale method)

1. Materials

1.1 Cell

- Cultured cells(Floating cells, Adhesive cells) >10⁹ cells

1.2 Reagents

- PBS(-) · 2-[4-(2-hydroxyethyl)-1-piperazinyl]ethanesulfonic acid (HEPES)
- Sodium hydroxide (NaOH) · Potassium chloride (KCl) · Magnesium chloride (MgCl₂)
- Sodium chloride (NaCl) · Calcium chloride (CaCl₂)
- Ethylenediamine tetraacetic acid (EDTA) · Glycerol (Glycerin) · Nonidet P-40 (NP-40)
- Dithiothreitol (DTT) · Phenyl fluoride methane sulfonyl (PMSF)

1.3 Apparatus

- High-speed cold centrifuge(HITACHI CF15RXII, BECKMAN COULTER Avanti HP301, etc.)
- Homogenizer(Kontes all glass Dounce homogenizer (B type pestle), etc.)
- Magnetic stirrer

2. Method

2.1 Preparation for reagent solutions

- 1) Buffer A: 10 mM HEPES-NaOH (pH7.9), 10 mM KCl, 1.5 mM MgCl₂, 0.5 mM DTT
- 2) Buffer B: 0.3 M HEPES-NaOH (pH7.9), 1.4 M KCl, 30 mM MgCl₂
- 3) Buffer C: 20 mM HEPES-NaOH (pH7.9), 0.42 M NaCl, 1.5 mM MgCl₂, 0.2 mM EDTA, 25%(v/v) glycerol, 0.5 mM DTT, 0.5 mM PMSF
- 4) Buffer D: 20 mM HEPES-NaOH (pH7.9), 100 mM KCl, 1 mM MgCl₂, 0.2 mM CaCl₂, 0.2 mM EDTA, 10%(v/v) glycerol, 1 mM DTT, 0.2 mM PMSF

2.2 Procedures

2.2.1 Harvest cells from cell culture media

For cells in suspension:

- 1) Collect cells in an appropriate centrifuge tube, and centrifuge at 500 ×g for 10 min at 4 °C. and remove the supernatant.
- 2) Suspend cells in PBS(-) and transfer to an appropriate centrifuge tube.
- 3) Centrifuge at 2,000 rpm for 5 min at 4 °C and remove the supernatant.

For adherent cells:

- 1) Wash cells with PBS(-), once or two times.
- 2) Add PBS(-), collect cells by scraper, and transfer to an appropriate centrifuge tube.
- 3) Centrifuge at 1,000 rpm for 5 min at 15 °C and remove the supernatant.

2.2.2 Crushing of whole cells

- 4) Wash with five packed cell volumes(5PCV) of PBS(-).
- 5) Centrifuge at 2,000 rpm for 5 min at 4 °C, and remove the supernatant.
- 6) Suspend in 5PCV of Buffer A, and stand on ice for 10 min.
- 7) Centrifuge at 2,000 rpm for 5 min at 4 °C, and remove the supernatant.
- 8) Suspend in 2PCV of Buffer A, and lyse by 10-20 strokes of a Kontes all glass Dounce homogenizer (B type pestle).
- 9) Check the homogenate microscopically for cell lysis.
- 10) Centrifuge at 600 ×g for 10 min at 4 °C, and transfer the supernatant (cytoplasmic extract) and the pellet (nuclei) to each flesh tube, and go to those methods.

2.2.3 Separation of cytoplasmic fraction and cell membrane fraction

- 11) Add 0.11PCV of Buffer B to the supernatant.
- 12) Centrifuge at 100,000 ×g for 60 min at 4 °C, and transfer the supernatant to a flesh tube.

Experiment Protocol 401

(Precipitation contains cell membrane fraction)

- 13) Dialyze against with 50 volumes of Buffer D at the collected supernatant at 4 °C.
- 14) Centrifuge at 100,000 ×g for 60 min at 4 °C, and transfer the supernatant to a flesh tube.
(Determine protein concentration by Bradford protein assay.)
- 15) Freeze with liquid nitrogen and store at -80 °C. [Cytoplasmic fraction]

2.2.4 Solubilization of cell membrane fraction

- 16) Add the surfactant to the precipitate of Step 12) so that the final concentration is as follows, and suspend by pipetting. Furthermore leave overnight at 4 °C.

Types of surfactants:

- a) 1 %-3 % Octylglucoside
 - b) 1 % CHAPS
 - c) 1 % NP-40
 - d) 1 % Tween20
 - e) 1 % Triton X-100
- 17) Octylglucoside is removed by dialysis. Other surfactants are diluted with Buffer D to a surfactant concentration of 0.1 %.
(Determine protein concentration by Bradford protein assay.)
 - 18) Freeze with liquid nitrogen and store at -80 °C. [Membrane fraction]

2.2.5 Preparation of nuclear extract

- 19) Centrifuge the pellet from Step 10) at 20,000 ×g for 20 min at 4 °C and remove the supernatant.
- 20) Suspend in 1PCV of Buffer C, and lyse by 10 strokes of a Kontes all glass Dounce homogenizer (B type pestle).
- 21) Stir gently with a magnetic stirring bar in beaker for 30 min at 4 °C.
- 22) Centrifuge at 20,000 ×g for 30 min at 4 °C, and transfer the supernatant to a flesh tube.
- 23) Dialyze against 50 volumes of Buffer D for 5h at 4 °C.
- 24) Centrifuge at 20,000 ×g for 30 min at 4 °C, and transfer the supernatant to a flesh tube.
(Determine protein concentration by Bradford protein assay.)
- 25) Freeze with liquid nitrogen and store at -80 °C. [Nuclear extract]

3. Supplements

- It is desirable that the amount of cells(PCV) before fractionation in Step 4) is about 5-10 ml.
- When using a homogenizer, pull it up slowly without foaming.
- If the detergent concentration is 1 % or higher, it may interfere with the binding of the ligand to the target protein during affinity purification. Reduce the surfactant concentration to 0.1 % before use.

References: J.D.Dignam, R.M.Lebovitz, and R.G.Roeder, Nucleic Acids Res. 11, 1475 (1983)