

Experiment Protocol 402

Preparation of cell extract (Small scale method)

1. Materials

1.1 Cell

- Cultured cells(Floating cells, Adhesive cells) $>10^7$ cells

1.2 Reagents

- PBS(-) · 2-[4-(2-hydroxyethyl)-1-piperazinyl]ethanesulfonic acid (HEPES)
- Sodium hydroxide (NaOH) · Potassium chloride (KCl) · Magnesium chloride ($MgCl_2$)
- Sodium chloride (NaCl) · Calcium chloride ($CaCl_2$)
- Ethylenediamine tetraacetic acid (EDTA) · Glycerol (Glycerin) · Nonidet P-40 (NP-40)
- Dithiothreitol (DTT) · Phenyl fluoride methane sulfonyl (PMSF)

1.3 Apparatus

- High-speed cold centrifuge(HITACHI CF15RXII, etc.)
- Vortex mixer

2. Method

2.1 Preparation for reagent solutions

- 1) Buffer A: 10 mM HEPES-NaOH (pH7.9), 10 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5 mM PMSF
- 2) Buffer B: 20 mM HEPES-NaOH (pH7.9), 0.42 M NaCl, 1 mM EDTA, 1 mM DTT, 0.5 mM PMSF
- 3) Buffer C: 20 mM HEPES-NaOH (pH7.9), 100 mM KCl, 1 mM $MgCl_2$, 0.2 mM $CaCl_2$, 0.2 mM EDTA, 10%(v/v) glycerol, 1 mM DTT, 0.2 mM PMSF

2.2 Procedures

2.2.1 Fractionation

For cells in suspension:

- 1) Collect cells in an appropriate centrifuge tube, and centrifuge at $500 \times g$ for 10 min at 4 °C. and remove the supernatant.
- 2) Suspend cells in PBS(-) and transfer to an appropriate centrifuge tube.
- 3) Centrifuge at 2,000 rpm for 5 min at 4 °C and remove the supernatant.

For adherent cells:

- 1) Wash cells with PBS(-), once or two times.
- 2) Add PBS(-), collect cells by scraper, and transfer to an appropriate centrifuge tube.
- 3) Centrifuge at 1,000 rpm for 5 min at 15 °C and remove the supernatant.

2.2.2 Lysis of whole cells

- 4) Wash with five packed cell volumes(5PCV) of PBS(-).
- 5) Centrifuge at $1,500 \times g$ for 5 min at 4 °C, and remove the supernatant.
- 6) Repeat the above Step 4) to 5) two more times. (Wash the cells with buffer three times in total.)
- 7) Suspend by pipetting in 2PCV of Buffer A.
- 8) Stand on ice for 15 min.
- 9) Add the surfactant to the precipitate so that the final concentration is as follows, and suspend by pipetting.
Types of surfactants:
 - a) 1 %-3 % Octylglucoside
 - b) 1 % CHAPS
 - c) 1 % NP-40
 - d) 1 % Tween20
 - e) 1 % Triton X-100
- 10) Vortex for 10 sec.
- 11) Centrifuge at $20,000 \times g$ for 5 min at 4 °C, and transfer the supernatant(cytoplasmic extract)

Experiment Protocol 402

and the pellet (nuclei) to each flesh tubes.

2.2.3 Separation of cytoplasmic fraction and cell membrane fraction

- 12) Octylglucoside is removed by dialysis. Other surfactants are diluted with Buffer C to a surfactant concentration of 0.1 %.
(Determine protein concentration by Bradford protein assay.)
- 13) Freeze with liquid nitrogen and store at -80 °C. [Mixture of cytoplasmic fraction and membrane fraction]

2.2.4 Preparation of nuclear extract

- 14) Suspend the pellet from Step 11) in 1PCV of Buffer B, and incubate for 15 min at 4 °C.
- 15) Centrifuge at 20,000 ×g for 5 min at 4 °C, and transfer the supernatant to a flesh tube.
- 16) Dialyze against 50 volumes of Buffer C at 4 °C.
- 17) Centrifuge at 20,000 ×g for 30 min at 4 °C, and transfer the supernatant to a flesh tube.
(Determine protein concentration by Bradford protein assay.)
- 18) Freeze with liquid nitrogen and store at -80 °C. [Nuclear extract]

3. Supplements

- If the detergent concentration is 1 % or higher, it may interfere with the binding of the ligand to the target protein during affinity purification. Reduce the surfactant concentration to 0.1 % before use.