## Immunoprecipitation (4)

### Introduction

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## Materials and method

### Materials

1. Protein G beads
2. *S. cerevisiae* cells expressing HA-Opt2
3. PBS(-)
4. Lysis buffer (PBS(-), 1 mM PMSF, 1 mM DTT, 1x protease inhibitor cocktail)
5. Wash buffer (10 mM HEPES-NaOH(pH7.9), 50 mM KCl, 0.2 mM EDTA, 10%(v/v) glycerol)
6. Anti-HA-tag antibody (TANA2; from MBL)
7. Anti-Pma1 antibody (yN-20; from Santa Cruz)
8. Anti-Dpm1 antibody (5C5A7; from Thermo Fisher)
9. Anti-Pho8 antibody (gift from Prof. Y. Ohsumi (Tokyo Institute of Technology))
10. Anti-Pep12 antibody (2C3G4; from Thermo Fisher)
11. HA peptide buffer (1 mg/mL HA peptide (from MBL) in PBS(-))
12. TBS-T buffer (TBS(-), 0.1% Tween 20)
13. Blocking buffer (5% skim milk in TBS-T buffer)
14. Anti-mouse antibody HRP conjugate (from GE Healthcare)
15. Anti-rabbit antibody HRP conjugate (from GE Healthcare)

### Method 1 (Lysate Preparation)

1. **Yeast culture**  
Culture yeast cells to log phase in YPD medium (1% yeast extract, 2% pepton, 2% glucose).
2. **Harvest and homogenize**  
Collect cells by centrifugation and wash with PBS(-). Resuspend cells with Lysis buffer, and lysis them by vigorously mixing with glass beads for 10 min at 4°C. Transfer the lysate to new tube without taking beads.
3. **Vesiculate membrane**  
Sonicate the lysate (10 sec x 3) to vesiculate the membrane. Centrifuge at 13,000 g for 15 min. Ultracentrifuge the supernatant at 100,000 g for 30 min. Remove the supernatant, resuspend the pellet with Lysis buffer, and sonicate (10 sec x 2) ("Vesiculated membrane").

### Method 2 (Immunoprecipitation)

1. **Beads preparation**  
Transfer 1 mg (50 ul) of beads to a tube. Wash beads with 200 ul ice-cold PBS(-) twice.
2. **Bind Antibody**  
Add antibody (50 ug in 200 ul PBS(-)) to beads. Mix for 30 min at room temperature. Wash beads bound to antibody with Wash buffer 3 times, and store at 4°C in 200 ul Wash buffer.
3. **Immunoprecipitation**  
Add 100 ul beads bound to the antibody to "Vesiculated membrane", and incubate overnight at 4°C with rotation.
4. **Wash**  
Separate beads magnetically from the suspension, and collect a part of supernatant for immunoblot analysis (Unbound fraction). Wash beads 3 times with PBS(-).  
Optional\*: Resuspend beads with HA peptide buffer and stand for 5 min at 4°C. Separate magnetically and collect supernatant.  
\*, In some cases, vesicles containing HA-tagged protein may be released from the beads during this short period of incubation. If so, omit this process.
5. **Elution**  
Resuspend beads with HA peptide buffer and incubate for 8 h at 4°C with rotation. Separate the beads magnetically and collect supernatant (Elution fraction).

### FG beads® information

Product name	Protein G beads
Product number	TAS8848N1173
Storage temperature	2-8°C
Storage buffer	10 mM HEPES(pH7.9), 50 mM KCl, 1 mM EDTA, 10% glycerol
Size of beads	190 nm ± 20 nm
IgG binding capacity	>100 ug mouse IgG /mg of beads

### Methods 3 (Western blotting)

1. Perform SDS-PAGE.
2. Transfer proteins from the gel to a PVDF membrane.
3. Block the membrane with Blocking buffer for 60 min at room temperature.
4. Dilute the each primary antibody with the Blocking buffer. (Anti-HA; 1/2,000, Anti-Pma1; 1/1,000, Anti-Dpm1; 1/500, Anti-Pho8; 1/1,000, Anti-Pep12; 1/1,000)
5. Incubate the membrane for 60 min at room temperature.
6. Wash the membrane with TBS-T buffer 3 times.
7. Diluted the Secondary antibody with TBS-T buffer to 1/7,500, and incubate the membrane for 60 min at room temperature.
8. Wash the membrane with TBS-T buffer 3 times.
9. Detect with a chemiluminescence substrate.

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