Immunoprecipitation (3)

Introduction

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
Materials and method

Materials

1. Protein G beads
2. *S.cerevisiae* lysate
3. Lysis buffer (20mM HEPES(pH8.0) , 9M Urea)
4. Protease inhibitor cocktail (from Nacalai Tesque)
5. PR-619 (Life Technologies)
6. Triton buffer (1% Triton-X100, 20mM Tris-HCl(pH7.4), 150mM NaCl, 2mM EDTA)
7. Anti-Ssa1 antibody (from Santa Cruz)
8. SDS sample buffer
9. Anti-mono- and polyubiquitinated conjugates antibody (FK2H) HRP-conjugate (from Enzo Life Sciences)
10. PBS-T buffer (137mM NaCl, 2.7mM KCl, 8.1mM Na₂HPO₄·12H₂O, 1.5mM KH₂PO₄, 0.1% Tween 20)
11. Blocking buffer (5% skim milk in PBS-T buffer)

FG beads® information

Product name	Protein G beads
Product number	TAS8848N1173
Storage temperature	2-8°C
Storage buffer	10mM HEPES(pH7.9), 50mM KCl, 1mM EDTA, 10%glycerol
Size of beads	190nm ± 20 nm
IgG binding capacity	>100ug mouse IgG /mg of beads

Method 1 (Lysate Preparation)

1. Yeast culture

Culture yeast for 3 days at 30°C in synthetic complete medium lacking leucine (SC-Leu).

2. Harvest and homogenize

Centrifuge at 3000 × g for 10min, wash twice in sterile water, and resuspend in lysis buffer with protease inhibitor cocktail and 50uM PR-619.
Vortex with glass beads and centrifuge to obtain lysate.

Method 2 (Immunoprecipitation)

1. Lysate dilution

Add 9 volume of Triton buffer to lysate.

2. Reaction

Incubate lysate containing 125ug of total protein and 0.25ug of Anti-Ssa1 antibody at room temperature for 60min on microtube mixer.

3. Collect antigen-antibody complex

Add 25ul of Protein G beads, incubate at room temperature for 60min on microtube mixer.

4. Wash

Separate magnetically and remove supernatant.
Wash beads with Triton buffer 5 times.
Separate magnetically and remove supernatant.

5. Elution

Add SDS sample buffer and resuspend beads.
Incubate for 10min at 70°C and remove the beads.

Methods 3 (Western blotting)

1. Perform SDS-PAGE.
2. Transfer the protein from the gel to a PVDF membrane using Trans-Blot® Turbo™.
3. Block the membrane with the Blocking buffer for 120min at room temperature.
4. Dilute the Anti-ubiquitine antibody with the Blocking buffer to 1/1000.
5. Incubate the membrane for 120min at room temperature.
6. Wash the membrane with PBS-T buffer 3 times.
7. Detect with a chemiluminescence substrate.

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